

Product Datasheet

COLLAGEN I/II/III/IV/V antibody (orb429972)

Description	Rabbit polyclonal antibody to COLLAGEN I/II/III/IV/V
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr
Immunogen	A mixture of human collagen I, II, III, IV and V from human placenta.
Target	COLLAGEN I/II/III/IV/V
Preservatives	0.1%Mannitol0.1%Dextran0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.