

Product Datasheet

COLLAGEN I/III antibody (orb429976)

Description	Rabbit polyclonal antibody to COLLAGEN I/III
Species/Host	Rabbit
Reactivity	Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr
Immunogen	Porcine collagen, types 1 and 3 from skin.
Target	COLLAGEN I/III
Preservatives	1% Mannitol1% Dextran. 0.1M Sodium Chloride0.01M Sodium Phosphate0.01M Sodium Borate
Form/Appearance	Purified IgG - lyophilized
Storage	-20°C after reconstitution
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.