

Product Datasheet

CORTICOTROPHIN RELEASING FACTOR antibody (orb429986)

Description	Mouse monoclonal antibody to CORTICOTROPHIN RELEASING FACTOR
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	Corticotrophin Releasing Factor conjugated to Bovine Serum Albumin.
Target	CORTICOTROPHIN RELEASING FACTOR
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified Ig - liquid
Concentration	1mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.