

Product Datasheet

CORTISOL antibody (orb429992)

Description	Sheep polyclonal antibody to CORTISOL
Species/Host	Sheep
Reactivity	Aves
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Cortisol-21-hemisuccinyl-KLH/BSA/Fetuin mixture.
Target	CORTISOL
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	5.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.