


Product Datasheet

CRYPTOSPORIDIUM antibody (orb429996)

Description	Mouse monoclonal antibody to CRYPTOSPORIDIUM
Species/Host	Mouse
Reactivity	Protozoa
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, ICC, IF
Immunogen	Purified cryptosporidium oocysts from bovine faeces.
Target	CRYPTOSPORIDIUM
Preservatives	PBS: 0.09% Sodium Azide (NaN ₃)
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.
Note	For research use only
Isotype	IgG3
Clonality	Monoclonal
Purity	Purified IgG prepared by affinity chromatography on Protein A from tissue culture supernatant.
Dilution Range	ELISA: 1/2000-8000, IF/ICC: 1/50-200
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
Cryptosporid...