

Product Datasheet

DIGOXIGENIN antibody (orb430018)

Description	Sheep polyclonal antibody to DIGOXIGENIN
Species/Host	Sheep
Reactivity	Chemical
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-Fr, WB
Immunogen	Digoxin-3-CM conjugated to keyhole limpet haemocyanin.
Target	DIGOXIGENIN
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	5 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.