

## Product Datasheet

### ENDOTHELIN 1 antibody (orb430036)

|                                |                                                                     |
|--------------------------------|---------------------------------------------------------------------|
| <b>Description</b>             | Mouse monoclonal antibody to<br>ENDOTHELIN 1                        |
| <b>Species/Host</b>            | Mouse                                                               |
| <b>Reactivity</b>              | Human                                                               |
| <b>Conjugation</b>             | Unconjugated                                                        |
| <b>Tested<br/>Applications</b> | IH-Fr, IRMA                                                         |
| <b>Immunogen</b>               | C-terminal fragment of human<br>endothelin 1.                       |
| <b>Target</b>                  | ENDOTHELIN 1                                                        |
| <b>Preservatives</b>           | 0.1% Sodium Azide (NaN <sub>3</sub> ). Phosphate<br>buffered saline |
| <b>Form/Appearance</b>         | Purified IgG - liquid                                               |
| <b>Concentration</b>           | 1.0 mg/ml                                                           |
| <b>Storage</b>                 | +4°C, -20°C if preferred                                            |
| <b>Note</b>                    | For research use only                                               |
| <b>Isotype</b>                 | IgG1                                                                |
| <b>Clonality</b>               | Monoclonal                                                          |
| <b>Purity</b>                  | Purified                                                            |
| <b>Expiration Date</b>         | 12 months from date of receipt.                                     |