

Product Datasheet

ESCHERICHIA COLI antibody (Biotin) (orb430052)

Description	Rabbit polyclonal antibody against ESCHERICHIA COLI conjugated to Biotin
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	A mixture of all antigenic serotypes.
Target	ESCHERICHIA COLI
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	4.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.