

Product Datasheet

ESCHERICHIA COLI J5 LPS antibody (orb430055)

Description	Mouse monoclonal antibody to ESCHERICHIA COLI J5 LPS
Species/Host	Mouse
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	E. coli J5 cells.
Target	ESCHERICHIA COLI J5 LPS
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.