

Product Datasheet

ESCHERICHIA COLI K88A antibody (orb430056)

Description	Mouse monoclonal antibody to ESCHERICHIA COLI K88A
Species/Host	Mouse
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF
Immunogen	E. coli.
Target	ESCHERICHIA COLI K88A
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.