

Product Datasheet

ESCHERICHIA COLI LABILE TOXIN A CHAIN antibody (orb430061)

Description	Mouse monoclonal antibody to ESCHERICHIA COLI LABILE TOXIN A CHAIN
Species/Host	Mouse
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Immunogen	Escherichia coli heat-labile enterotoxin.
Target	ESCHERICHIA COLI LABILE TOXIN A CHAIN
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.