

Product Datasheet

FELINE IMMUNODEFICIENCY VIRUS TRANSMEMBRANE PROTEIN gp41 antibody (orb430070)

Description	Mouse monoclonal antibody to FELINE IMMUNODEFICIENCY VIRUS TRANSMEMBRANE PROTEIN gp41
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Native, FIV extracts.
Target	FELINE IMMUNODEFICIENCY VIRUS TRANSMEMBRANE PROTEIN gp41
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.