

Product Datasheet

GLUCOSE TRANSPORTER 4 antibody (FITC) (orb430112)

Description	Mouse monoclonal antibody against GLUCOSE TRANSPORTER 4 conjugated to FITC
Species/Host	Mouse
Reactivity	Rat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	Partially purified vesicles containing insulin-responsive glucose transporter 4.
Target	GLUCOSE TRANSPORTER 4
Preservatives	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.