

Product Datasheet

GLUCOSE TRANSPORTER 5 antibody (orb430115)

Description	Rabbit polyclonal antibody to GLUCOSE TRANSPORTER 5
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	A synthetic peptide (ELKELPPVTSEQ) corresponding to the carboxy terminus of the human GLUT 5. To accomodate KLH conjugation a cysteine was added to the N-terminal.
Target	GLUCOSE TRANSPORTER 5
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	5.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.