

## Product Datasheet

### GLUTAMATE DEHYDROGENASE antibody (orb430119)

|                            |                                                                     |
|----------------------------|---------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to<br>GLUTAMATE DEHYDROGENASE            |
| <b>Species/Host</b>        | Rabbit                                                              |
| <b>Reactivity</b>          | Bovine                                                              |
| <b>Conjugation</b>         | Unconjugated                                                        |
| <b>Tested Applications</b> | ELISA, WB                                                           |
| <b>Immunogen</b>           | Bovine glutamate dehydrogenase<br>from liver.                       |
| <b>Target</b>              | GLUTAMATE DEHYDROGENASE                                             |
| <b>Preservatives</b>       | 0.1% Sodium Azide (NaN <sub>3</sub> ).<br>Phosphate buffered saline |
| <b>Form/Appearance</b>     | Purified Ig - liquid                                                |
| <b>Concentration</b>       | 10.0 mg/ml                                                          |
| <b>Storage</b>             | +4°C, -20°C if preferred                                            |
| <b>Note</b>                | For research use only                                               |
| <b>Isotype</b>             | Polyclonal IgG                                                      |
| <b>Clonality</b>           | Polyclonal                                                          |
| <b>Purity</b>              | Purified                                                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                     |