

Product Datasheet

GREEN FLUORESCENT PROTEIN antibody (orb430127)

Description	Sheep polyclonal antibody to GREEN FLUORESCENT PROTEIN
Species/Host	Sheep
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF
Immunogen	Green fluorescent protein from Aequorea victoria.
Target	GREEN FLUORESCENT PROTEIN
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	5.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.