

Product Datasheet

HUMAN IMMUNODEFICIENCY VIRUS 1 p24 antibody (orb430164)

Description	Mouse monoclonal antibody to HUMAN IMMUNODEFICIENCY VIRUS 1 p24
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, WB
Immunogen	HIV p24.
Target	HUMAN IMMUNODEFICIENCY VIRUS 1 p24
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.