

Product Datasheet

HUMAN IMMUNODEFICIENCY VIRUS 1 p24 antibody (orb430165)

Description	Goat polyclonal antibody to HUMAN IMMUNODEFICIENCY VIRUS 1 p24
Species/Host	Goat
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Purified native p24, from strain IIIB.
Target	HUMAN IMMUNODEFICIENCY VIRUS 1 p24
Preservatives	0.1% Sodium azide (NaN ₃).. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	4.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.