

Product Datasheet

LUTEINIZING HORMONE RELEASING HORMONE RECEPTOR antibody (orb430227)

Description

Mouse monoclonal antibody to

RELEASING HORMONE
RECEPTOR

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FC, IH-Fr, WB

Applications
Immunogen

Synthetic peptide corresponding to the amino acids 1 - 29 of the human luteinizing hormone releasing hormone receptor extracellular domain, conjugated to Bovine Serum Albumin.

Target

LUTEINIZING HORMONE
RELEASING HORMONE
RECEPTOR

Preservatives

0.09% Sodium Azide (NaN₃).
Phosphate buffered saline

Form/Appearance

Purified IgG - liquid

Concentration

1.0 mg/ml

Storage

+4°C, -20°C if preferred

Note

For research use only

Isotype

IgG1

Clonality

Monoclonal

Purity

Purified

Expiration Date

12 months from date of receipt.