

Product Datasheet

MELANOCYTE CELL SURFACE ANTIGEN antibody (orb430236)

Description	Mouse monoclonal antibody to MELANOCYTE CELL SURFACE ANTIGEN
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr
Immunogen	Human KW-mel(2.1) melanoma cells.
Target	MELANOCYTE CELL SURFACE ANTIGEN
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.