

Product Datasheet

ALPHA 1 MICROGLOBULIN antibody (orb430242)

Description	Sheep polyclonal antibody to ALPHA 1 MICROGLOBULIN
Species/Host	Sheep
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ID
Immunogen	Highly pure human alpha-1-microglobulin, prepared from the urine of patients with chronic renal proteinuria, molecular weight 25-30 kDa.
Target	ALPHA 1 MICROGLOBULIN
Preservatives	0.09% Sodium Azide (NaN ₃)>. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	5mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.