

Product Datasheet

BETA 2 MICROGLOBULIN antibody (orb430246)

Description	Mouse monoclonal antibody to BETA 2 MICROGLOBULIN
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Beta-2 microglobulin from human urine.
Target	BETA 2 MICROGLOBULIN
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.