

Product Datasheet

NEOMYCIN antibody (orb430275)

Description	Rabbit polyclonal antibody to NEOMYCIN
Species/Host	Rabbit
Reactivity	Chemical
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Neomycin conjugated to Keyhole limpet Haemocyanin.
Target	NEOMYCIN
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	5.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.