

Product Datasheet

17-BETA-ESTRADIOL antibody (orb430293)

Description	Mouse monoclonal antibody to 17-BETA-ESTRADIOL
Species/Host	Mouse
Reactivity	Aves
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	17-beta estradiol conjugated to bovine serum albumin at position 6.
Target	17-BETA-ESTRADIOL
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.