

Product Datasheet

PEPSINOGEN I antibody (orb430315)

Description	Mouse monoclonal antibody to PEPSINOGEN I
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, RIA
Immunogen	Purified human pepsinogen 1.
Target	PEPSINOGEN I
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.