

Product Datasheet

PEPTIDOGLYCAN antibody (orb430318)

Description	Mouse monoclonal antibody to PEPTIDOGLYCAN
Species/Host	Mouse
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-Fr, IHC-P
Immunogen	This antibody was raised against insoluble peptidoglycan obtained by TCA-heat and ethanol extraction of Streptococcus mutans BHT cells.
Target	PEPTIDOGLYCAN
Preservatives	None present.
Form/Appearance	Ascitic Fluid - raw
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Ascites
Expiration Date	12 months from date of receipt.