

Product Datasheet

PICHIA PASTORIS antibody (orb430323)

Description	Rabbit polyclonal antibody to PICHIA PASTORIS
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, WB
Immunogen	P. pastoris whole cells.
Target	PICHIA PASTORIS
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	4.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.