

Product Datasheet

PITUITARY ADENYLATE CYCLASE ACTIVATING POLYPEPTIDE antibody (orb430324)

Description	Sheep polyclonal antibody to PITUITARY ADENYLATE CYCLASE ACTIVATING POLYPEPTIDE
Species/Host	Sheep
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Proprietary, 14 amino acid peptide sequence of PACAP-38, conjugated to KLH.
Target	PITUITARY ADENYLATE CYCLASE ACTIVATING POLYPEPTIDE
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	5.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.