

Product Datasheet

PSEUDOMONAS AERUGINOSA antibody (orb430367)

Description	Mouse monoclonal antibody to PSEUDOMONAS AERUGINOSA
Species/Host	Mouse
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified outer membrane protein of Pseudomonas aeruginosa.
Target	PSEUDOMONAS AERUGINOSA
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.