

Product Datasheet

RESPIRATORY SYNCYTIAL VIRUS FUSION PROTEIN antibody (orb430378)

Description	Mouse monoclonal antibody to RESPIRATORY SYNCYTIAL VIRUS FUSION PROTEIN
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Purified Respiratory Syncytial Virus (RSV).
Target	RESPIRATORY SYNCYTIAL VIRUS FUSION PROTEIN
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.