

Product Datasheet

SALMONELLA ENTERITIDIS LPS antibody (orb430390)

Description	Mouse monoclonal antibody to SALMONELLA ENTERITIDIS LPS
Species/Host	Mouse
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Heat inactivated cells from Salmonella (LPS groups A, B & D).
Target	SALMONELLA ENTERITIDIS LPS
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	GROUP D 0-12 ANTIGEN
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.