

Product Datasheet

SEROTONIN antibody (orb430395)

Description	Rat monoclonal antibody to SEROTONIN
Species/Host	Rat
Reactivity	Aves
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr
Immunogen	Serotonin conjugated to bovine serum albumin.
Target	SEROTONIN
Preservatives	0.09% Sodium Azide (NaN ₃)
Form/Appearance	Tissue Culture Supernatant - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2
Clonality	Monoclonal
Purity	S/N
Expiration Date	12 months from date of receipt.