

Product Datasheet

SEROTONIN TRANSPORTER antibody (orb430396)

Description	Mouse monoclonal antibody to SEROTONIN TRANSPORTER
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IH-Fr
Immunogen	Rat serotonin transporter, N-terminus/GST fusion protein (aa 1-85).
Target	SEROTONIN TRANSPORTER
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(N-TERMINAL)
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.