

Product Datasheet

STREPTOCOCCUS B antibody (orb430413)

Description	Mouse monoclonal antibody to STREPTOCOCCUS B
Species/Host	Mouse
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Streptococcus B organisms types Ia, Ib, Ic, II and III.
Target	STREPTOCOCCUS B
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.8 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.