

Product Datasheet

STREPTOCOCCUS PNEUMONIAE antibody (orb430416)

Description	Mouse monoclonal antibody to STREPTOCOCCUS PNEUMONIAE
Species/Host	Mouse
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF
Immunogen	Streptococcus pneumoniae whole cells.
Target	STREPTOCOCCUS PNEUMONIAE
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.