

Product Datasheet

SUBSTANCE P antibody (orb430421)

Description	Rat monoclonal antibody to SUBSTANCE P
Species/Host	Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, IHC-P
Immunogen	Substance P conjugated to Bovine Serum Albumin.
Target	SUBSTANCE P
Preservatives	0.09% Sodium Azide (NaN ₃)
Form/Appearance	Tissue Culture Supernatant - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	S/N
Expiration Date	12 months from date of receipt.