

Product Datasheet

THYROID STIMULATING HORMONE antibody (orb430444)

Description	Mouse monoclonal antibody to THYROID STIMULATING HORMONE
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Human TSH beta subunit.
Target	THYROID STIMULATING HORMONE
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	BETA CHAIN
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.