

Product Datasheet

UBIQUITIN antibody (orb430475)

Description	Sheep polyclonal antibody to UBIQUITIN
Species/Host	Sheep
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Native ubiquitin, purified from bovine erythrocytes to a purity $\geq 90\%$ by SDS-PAGE.
Target	UBIQUITIN
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	5.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.