

Product Datasheet

ARG-VASOPRESSIN antibody (orb430482)

Description	Rabbit polyclonal antibody to ARG-VASOPRESSIN
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	Synthetic arginine vasopressin conjugated to bovine thyroglobulin
Target	ARG-VASOPRESSIN
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	5.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.