

Product Datasheet

YELLOW FEVER VIRUS antibody (orb430491)

Description	Mouse monoclonal antibody to YELLOW FEVER VIRUS
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	FA, ICC, IF
Immunogen	17D strain of yellow fever virus
Target	YELLOW FEVER VIRUS
Preservatives	None Present
Form/Appearance	Ascites - Liquid
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Ascites
Expiration Date	12 months from date of receipt.