

Product Datasheet

FGF 9 antibody (orb430502)

Description	Rabbit polyclonal antibody to FGF 9
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FA, WB
Immunogen	Recombinant murine FGF-9.
Target	FGF 9
Preservatives	None present.. Phosphate buffered saline
Form/Appearance	Purified IgG - lyophilized
Concentration	1.0 mg/ml
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.