

Product Datasheet

C-10 antibody (orb430505)

Description	Rabbit polyclonal antibody to C-10
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant mouse C-10.
Target	C-10
Preservatives	None Present. Phosphate buffered saline
Form/Appearance	Purified IgG - lyophilized
Concentration	1.0 mg/ml
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.