

Product Datasheet

G-CSF antibody (orb430509)

Description	Rabbit polyclonal antibody to G-CSF
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FA, WB
Immunogen	Recombinant mouse G-CSF.
Target	G-CSF
Preservatives	None Present. Phosphate buffered saline
Form/Appearance	Purified IgG - lyophilized
Concentration	1.0 mg/ml
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.