

Product Datasheet

G-CSF antibody (Biotin) (orb430510)

Description	Rabbit polyclonal antibody against G-CSF conjugated to Biotin
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Recombinant mouse G-CSF.
Target	G-CSF
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - lyophilized
Concentration	50 ug/ml
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.