

Product Datasheet

IP-10 antibody (Biotin) (orb430527)

Description	Rabbit polyclonal antibody against IP-10 conjugated to Biotin
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Recombinant mouse IP-10.
Target	IP-10
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - lyophilized
Concentration	0.1 mg/ml
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.