

Product Datasheet

MIP-1 GAMMA antibody (orb430538)

Description	Rabbit polyclonal antibody to MIP-1 GAMMA
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant mouse MIP-1 gamma.
Target	MIP-1 GAMMA
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - lyophilized
Concentration	1.0 mg/ml
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.