

Product Datasheet

VEGF antibody (orb430545)

Description	Rabbit polyclonal antibody to VEGF
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FA, ICC, IF, IHC-P, WB
Target	VEGF
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - lyophilized
Concentration	1.0 mg/ml
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.