

Product Datasheet

SF-20 antibody (Biotin) (orb430550)

Description	Rabbit polyclonal antibody against SF-20 conjugated to Biotin
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Recombinant mouse SF-20.
Target	SF-20
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - lyophilized
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.