

Product Datasheet

CDw198 antibody (orb430563)

Description	Goat polyclonal antibody to CDw198
Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthetic peptide PNVTMTDYYPDFFTAPCDAEFLLRGSM corresponding to amino acids 7-33 of mouse CDw198.
Target	CDw198
Preservatives	0.1% Sodium Azide (NaN ₃)0.1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.