

Product Datasheet

CD197 antibody (orb430565)

Description	Goat polyclonal antibody to CD197
Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, IHC-P, WB
Immunogen	Synthetic peptide DPGKPRKNVLVVALLVIFQVC corresponding to amino acids 2-22 of mouse CD197.
Target	CD197
Preservatives	0.1% Sodium Azide, 0.1% Bovine Serum Albumin.. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.